



## Dynamic PIP

### Feature Benefit

- ROM: PIP Extension: 65° Flexion to 25° Hyperextension  
PIP Flexion: 10° Flexion to 100° Flexion
- Bilateral tensioning system
- Tension adjustable dynamic springs for ideal load
- Length adjustable for optimal lever arm
- Circumference adjustable to accommodate edema
- Adjustable length to accommodate all fingers
- Available in extension or flexion

### Indications

- Crush Injury
- Tendon & Ligament Repair
- MUA
- Dupuytren's Contracture
- Fracture & Dislocation
- ORIF
- Burn
- Joint Arthroplasty
- Other Orthopedic and Neurological Conditions

