

Stat-A-Dyne™ WHFO

Wrist Hand Finger Orthosis





Stat-A-Dyne™ WHFO

(Wrist Hand Finger Orthosis)

Core Technologies

- Three modalities for tissue elongation (*static progressive, dynamic, & dual function*)
- Bi-directional hinge allows both flexion and extension in a single device
- Circular hinge for terminal range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment

Feature Benefit

- ROM: Extension 0° - 90° | Flexion 0° - 90° | Full Composite Fist
- Combined motion therapy for wrist and fingers
- Long Extrinsic Extensor/Flexor Stretch and Tenodesis Stretch
- Malleable finger attachments for individualized finger ROM
- Off the shelf or custom cuffs for maximum soft tissue capture
- Small, medium, and large gloves

Indications

- Fracture & Dislocation
- Sub-acute Tenolysis
- Crush Injury
- Tendon & Ligament Repair
- Burn
- Degloving
- CRPS
- Finger and/or Wrist Contracture
- Closed Manipulation
- Other Orthopedic and Neurological Conditions

Contact your local Territory Sales Consultant to order • www.lantzmedical.com • 866.236.8889