




## Vector 1 Hand CPM





## Vector 1 Hand CPM

### Feature Benefit

- Force gradable to accommodate the rehabilitation continuum
- Donning and doffing ease provided by unique glove system
- Range of motion beyond anatomical norm to overcome soft tissue and mechanical give
- ROM: -21° Hyperextension / 340° Flexion (full composite fist)
- Malleable finger springs to accommodate various ROM considerations
- Moveable finger plates allow joint isolation and safe treatment
- Fully programmable hand controller with compliance meter
- Capable of thumb adduction/abduction

### Indications

- Fracture & Dislocation
- Tenolysis
- Capsulotomy
- Crush Injury
- Tendon & Ligament Repair
- Burn
- Degloving
- CRPS
- Closed Manipulation
- Other Orthopedic and Neurological Conditions