

Stat-A-Dyne™ Knee Slider

Core Technologies

- Three modalities for tissue elongation (*static progressive, dynamic, & dual function*)
- Bi-directional hinge for flexion and extension
- Circular hinge for terminal end range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment
- Slider technology prevents joint compression & point loading

Feature Benefit

- ROM: -30° Hyperextension to 150° Flexion
- Bivalved or segmented Kydex cuffs maximize soft tissue capture and comfort
- Adjustable valgus/varus to accommodate male/female alignment
- Tibial tuberosity counterforce bracket provides 3 point lever arm maintaining joint alignment
- Long lever arm facilitates surface area distribution
- Foot/ankle cuff maintains anatomical alignment and provides 2 joint stretch

Indications

- Joint Arthroplasty
- Meniscus Repair
- Tibial Plateau Fracture
- Tendon & Ligament Repair
- ORIF
- Burn
- MUA
- Meniscectomy
- Other Orthopedic and Neurological Conditions

