



Stat-A-Dyne™ Pronation/Supination

Core Technologies

- Three modalities for tissue elongation (*static progressive, dynamic, & dual function*)
- Bi-directional hinge allows both pronation & supination in a single device
- Circular hinge for terminal end range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment

Feature Benefit

- Infinite ROM eliminates bottom out
- Bivalved or segmented Kydex cuffs maximize soft tissue capture and comfort
- Radial-based forearm cuff assures anatomically correct rotation
- Light-weight and low profile

Indications

- Fracture & Dislocation
- ORIF
- Crush Injury
- Radial Head Fracture
- Joint Arthroplasty
- Tendon & Ligament Repair
- Burn
- Closed Manipulation
- Other Orthopedic and Neurological Conditions

