



THUASNE



Vector 2
Hand CPM

Vector 2 Hand CPM

Recommended for: Fracture & Dislocation / Tenolysis / Capsulotomy / Crush Injury / Tendon & Ligament Repair / Burn / Degloving / CRPS / Closed Manipulation / Other Orthopedic and Neurological Conditions



Intuitive touch screen interface allows patient and provider to quickly adjust the device to address clinical presentation.

Feature Benefit

- Force gradable to accommodate the rehabilitation continuum
- Donning and doffing ease provided by unique glove system
- Range of motion beyond anatomical norm to overcome soft tissue and mechanical give
- ROM: -21° Hyperextension / 339° Flexion (280° full composite fist)
- Malleable finger springs to accommodate various ROM considerations
- Moveable finger plates allow joint isolation and safe treatment
- Fully programmable hand controller with compliance meter
- Capable of thumb adduction/abduction



Hyperextension capable to address palmer contractures



Detachable finger shoes may be placed on specific phalanges



Active tendon excursion without removing the product increases compliance and outcomes



Actuator bar removes for device servicing

Lantz Medical is Thuasne USA

*New brand, same great products and services,
now with more options than ever before!*

Lantz Medical is now Thuasne (Tu-An) USA, the American subsidiary of The Thuasne Group. Founded in 1847, The Thuasne Group is a medical device manufacturer based in Saint-Etienne, France with subsidiaries across the globe. Lantz Medical joins American medical device manufacturers Townsend Design and Quinn Medical under the Thuasne USA banner.

What changes for me?

With the new brand comes more products than ever before. As Thuasne USA, you'll be able to prescribe the same rehabilitation splinting you know and love while gaining access to a wide variety of orthotic bracing.

Thuasne USA Orthotics:

- Post-Operative – A wide variety of options including knee, back, elbow and more!
- Knee Bracing – Ligament/functional and osteoarthritis
- Back Bracing – LSO and TLSO, modular and stepdown styling that can be molded and modified for your patient
- Ankle Foot Orthoses – Composite AFO's in a variety of shell configurations and sizes with custom and off-the-shelf options



It is our pleasure to work with you to provide innovative and effective solutions to patients around the world. Your business is important to us and we are here with you every step of the way. If you have questions regarding the rebrand or need any assistance during this time of transition, please speak with your Thuasne USA representative or call 800.432.3466



Thuasne USA

4615 Shepard Street,
Bakersfield, CA 93313
800.432.3466

ThuasneUSA.com